


## AG-SILI TUBE BLU



Consists of a tubular fabric support in stretch mesh completely covered by polymer gel containing mineral oil and Silver ions; fits perfectly the toe in its shape, wraps it completely, softens irritated skin thanks to the soothing effect of mineral oil, is ideal for toe protection.

Is particularly suitable to protect callus, hyperkeratosis, to prevent skin sore due to rubbing or excessive pressure from the shoes; can be cut to desired size with common scissors.

Indications: toe callosity, corns, blisters, hammertoe, rubbing frices and excessive pressure.

Mis.: S, L

Cod.: UK009BF